

**IMPORTANT SAFETY INFORMATION**  
**EU Commercial Products**

## MIVI Q™ Aspiration Catheters

- Q3-36163-E, Q4-43150-E, Q5-57145-E, Q6-69145-E
- See IFU packaged with product for complete instructions on device usage (ref. IFU 101546.E)

**Intended Use/Indications for Use**

The Q Aspiration Catheter is indicated for the removal of fresh, soft emboli and thrombi in the peripheral and neurovascular systems. It may also be used as a diagnostic angiographic catheter.

**Contraindications**

- Vessels smaller than the catheter's distal outer diameter.

**Warnings**

- The Q Aspiration Catheter is intended to be used by physicians skilled in the field(s) of Interventional Radiology, Interventional Neuroradiology, Neurosurgery and/or Neurology who are skilled in endovascular interventional procedures for acute ischemic stroke or acute limb ischemia.
- Do not use the Q Aspiration Catheter with a guide catheter/sheath < 90 cm in length. This may result in the catheter body extending past the guide/sheath.
- Do not use the device in patients who are allergic to polyether block amides (PEBAX), polytetrafluoroethylene (PTFE) or polyether based thermoplastic polyurethanes (TPU).
- Never advance or retract the catheter against resistance without careful assessment of cause using fluoroscopy. If cause cannot be determined, withdraw the catheter.
- Do not torque the control wire of the Q Aspiration Catheter. Torqueing, movement against resistance or forced insertion of the catheter may result in device damage, vessel damage and/or tip breakage which may result in injury to the patient.
- If flow through catheter becomes restricted, do not attempt to clear catheter lumen by infusion. Doing so may cause catheter damage or patient injury. Remove or replace catheter.
- Use caution when injecting contrast agents into the distal cerebral neurovasculature through the guide catheter/sheath with a Q Aspiration Catheter in place. The potential for an increase in flow through the Q Aspiration Catheter may result in injury to the patient.

**Precautions**

- Store in a dark, dry place.
- Use prior to the "Use By" date.
- Do not expose device to solvents.

**IMPORTANT SAFETY INFORMATION**  
**EU Commercial Products**

- 
- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or reuse may result in impaired structural integrity and/or function.
  - Do not use open or damaged packages.
  - Upon removal from package, inspect device to ensure it is not damaged. Do not use kinked or damaged devices.
  - Confirm the compatibility of the Q Aspiration Catheter with any interventional devices to be used in conjunction with the catheter (e.g. guide catheter, microcatheter, guidewire).
  - Exercise care when handling the catheter before and during a procedure to reduce the possibility of damage.
  - Upon insertion and removal of the catheter, ensure the hemostasis valve is large enough to accommodate smooth movement.
  - Manipulate the catheter under fluoroscopic visualization. Do not attempt to move the catheter without observing the resultant tip response.
  - Avoid angular strain on the control wire during preparation and insertion.
  - Ensure that the proximal radiopaque marker does not extend past the distal marker band of the guide catheter.
  - When performing aspiration, ensure that the stopcock on the rotating hemostasis Y valve is open for only the minimum time needed to remove thrombus. Excess aspiration or failure to close the stopcock when aspiration is complete is not recommended.
  - Excessive force on the Q Aspiration Catheter is not recommended. If repositioning of the catheter is necessary, repositioning should be performed using appropriate techniques.
  - Do not exceed 300 kPa during contrast injection.
  - Do not use automated high-pressure contrast injection equipment with the Q Aspiration Catheter because it may damage the device.
  - If an intraluminal device becomes lodged in catheter, or if the catheter becomes severely kinked, withdraw the entire system (intraluminal device, catheter and guiding catheter).

**Adverse Events**

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with possible complications. Potential complications include, but are not limited to, the following:

- access site complications (hemorrhage, hematoma, fistula, etc.)
- air embolism
- allergic reaction (contrast, materials, device)
- death



6545 City West Parkway  
Eden Prairie, MN 55344 USA  
+1-952-944-3834

[www.mivineuro.com](http://www.mivineuro.com)

[customerservice@mivineuro.com](mailto:customerservice@mivineuro.com)

**IMPORTANT SAFETY INFORMATION**  
**EU Commercial Products**

---

- device malfunction
- emboli/distal embolization
- inability to remove all thrombus
- infection
- intracranial hemorrhage
- ischemia
- kidney damage
- neurological deficits or new stroke
- vasospasm
- vessel trauma (dissection or perforation)

**IMPORTANT SAFETY INFORMATION**  
**EU Commercial Products**

---

## MIVI Super 90™ 8F Guide Catheters

- MIA-9080S-E, MIA-9095S-E
- See IFU packaged with product for complete instructions on device usage (ref. IFU 101214.E).

### Intended Use/Indications for Use

The Super 90 8F Guide Catheter is indicated for use in facilitating the insertion and guidance of catheters into a selected blood vessel in the peripheral and neuro vascular system. It may also be used as a diagnostic angiographic catheter.

### Contraindications

None known

### Warnings

- The Super 90 8F Guide Catheter should only be used by physicians who are skilled in interventional procedures in the peripheral and neuro vascular systems.
- Do not use the device in patients who are allergic to polyether block amides (PEBAX/ VESTAMID), polytetrafluoroethylene (PTFE), polycarbonate or polyolefins.
- Never advance or retract the catheter against resistance without careful assessment of the cause using fluoroscopy. If a cause cannot be determined, withdraw catheter. Excessive torqueing, movement against resistance or forced insertion may result in catheter damage, vessel damage and/or tip breakage which may result in injury to the patient.
- If flow through catheter becomes restricted, do not attempt to clear catheter lumen by infusion. Doing so may cause catheter damage or patient injury. Remove or replace catheter.

### Precautions

- Store in a dark, dry place.
- Use prior to the "Use By" date.
- Do not expose device to solvents.
- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or reuse may result in impaired structural integrity and/or function.
- Confirm the compatibility of the Super 90 8F Guide Catheter with any interventional devices to be used in conjunction with the catheter.
- Do not use open or damaged packages.
- Upon removal from package, inspect device to ensure it is not damaged. Do not use kinked or damaged devices.



6545 City West Parkway  
Eden Prairie, MN 55344 USA  
+1-952-944-3834

[www.mivineuro.com](http://www.mivineuro.com)

[customerservice@mivineuro.com](mailto:customerservice@mivineuro.com)

## **IMPORTANT SAFETY INFORMATION**

### **EU Commercial Products**

---

- Exercise care when handling the catheter before and during a procedure to reduce the possibility of damage.
- Once hydrated, do not allow the device to dry out.
- Manipulate the catheter under fluoroscopic visualization. Do not attempt to move the catheter without observing the resultant tip response.
- Do not use automated high-pressure contrast injection equipment with the catheter, it may damage the device.
- Excessive force on the Super 90 8F Guide Catheter is not recommended. If repositioning of the catheter is necessary, repositioning should be performed using appropriate techniques.
- If intraluminal device becomes lodged in catheter, or if the catheter becomes severely kinked, withdraw the entire system (intraluminal device, catheter and guiding catheter).

#### **Potential Complications**

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with possible complications. Potential complications include, but are not limited to, the following:

- Access site complications (hemorrhage, hematoma, AV fistula, etc.)
- Allergic reaction (to contrast, materials or device)
- Death
- Embolism (air, thrombus or device)
- Infection
- Intracranial hemorrhage
- Ischemia
- Kidney damage
- Neurological deficits including stroke
- Vasospasm
- Vessel trauma (dissection or perforation)



6545 City West Parkway  
Eden Prairie, MN 55344 USA  
+1-952-944-3834

[www.mivineuro.com](http://www.mivineuro.com)

[customerservice@mivineuro.com](mailto:customerservice@mivineuro.com)

**IMPORTANT SAFETY INFORMATION**  
**EU Commercial Products**

---

## MIVI High Flow Tubing

- **HFT 110-E**
- See IFU packaged with product for complete instructions on device usage (ref. IFU 100813.F)

### **Intended Use/Indications for Use**

MIVI's High Flow Tubing is intended to be used in conjunction with commercially available aspiration pumps and aspiration catheters to aspirate, remove, or sample body fluids.

### **Contraindications**

None known

### **Warnings**

None

### **Precautions**

- Use prior to the "Use By" date.
- The device is intended for single use only. Do not resterilize or reuse.